





Radial Transmitters - Enthesopathies/Trigger Points





- Ø 6 mm
- Shock wave acupuncture
- Penetration depth: 0 20 mm
- Intensity:

F15





- Ø 15 mm
- Superficial pain regions, muscles of mastication, cervical spine
- Penetration depth: 0 30 mm
- Intensity:

T10 »Finger«

- Ø 10 mm
- Trigger points, deep-lying enthesopathy
- Penetration depth: 0 40 mm
- Intensity:

C15



»Basic Energy«



000000000

- Ø 15 mm
- Any type of enthesopathy
- Penetration depth: 0 35 mm
- Intensity:



»Deep Impact®«

DI15

- Ø 15 mm Deep-lying pathologies, chronic conditions, high-energy application
- Penetration depth: 0 50 mm
- Intensity:



•••••

000000000

Ø 15 mm

»CERAma-x®«

- Any type of enthesopathy
- Penetration depth: 0 30 mm
- Intensity:



000000000

000000000

D-ACTOR® Transmitters – Myofascial Treatment

D20-S »Basic Oscillator«



000000000

•••••

- Ø 20 mm
- Muscle tissue and connective tissue
- Penetration depth: 0 50 mm
- Intensity:



D20-T

Ø 20 mm For professional athletes. high-energy application

»Golden Oscillator«

- Penetration depth: 0 50 mm
- Intensity:



000000000

- D35-S
- »Large Oscillator«





- Muscle tissue and connective tissue. large muscles
- Penetration depth: 0 60 mm
- Intensity:
- **-**000000000

Special Transmitters – Enthesopathies/Trigger Points

Ro40





B15



- - **ATLAS** »Soft Tip« (after Di Maio/Swart)



- Ø 15 mm
- Any type of enthesopathy ■ Penetration depth: 0 – 35 mm
- Intensity:
- -----
- »White Tip«



- Any type of enthesopathy
- Intensity:
- Penetration depth: 0 30 mm

Ø 15 mm

- Any type of enthesopathy
- Penetration depth: 0 10 mm
- Intensity:
- •00000000



















SPINE-ACTOR® Transmitters – Treatment of the Paravertebral Muscles







The SPINE-ACTOR® transmitters have been developed to treat the paravertebral muscles. They come in three sizes allowing treatment in all sections of the spinal column muscles (lumbar, thoracic and cervical).

