

# CARDIOVIT AT-102 G2

Easy, smart and state-of-the-art



**SCHILLER**  
The Art of Diagnostics

# CARDIOVIT AT-102 G2

## The new generation ECG: CARDIOVIT AT-102 G2

- ⚡ Best signal quality and diagnostic algorithms for adult and paediatric patients
- ⚡ Fast bidirectional Wi-Fi communication with HIS/EMR
- ⚡ Spirometry and basic exercise ECG optionally available



SpiroScout SP plus, for spirometry tests with the CARDIOVIT AT-102 G2

## STATE-OF-THE-ART TECHNOLOGY AND CONNECTIVITY

The CARDIOVIT AT-102 G2 features an 8" colour display with touch function keys. Its sealed alphanumeric keyboard ensures easy cleaning, making it highly suitable for the daily clinical routine. Smart battery technology guarantees more than 8 hours of ECG recording.

Bidirectional communication allows for easy data access and fast transmission of all reports to EMR/HIS systems in PDF, DICOM or HL7 format. Standard Wi-Fi and LAN connectivity, combined with ECG review on the display and a large memory, supports paperless workflows and cost saving. Customizable printouts are available on the internal high-resolution A4 thermal printer.

## HIGH EFFICIENCY IN ANY WORKFLOW

ECG acquisition has never been easier and faster: Enter or download patient data with PDQ or worklist, acquire ECG, save, print and transmit data within twenty seconds using only one navigation key.




## EASY 1-2-3 STEPS

- 1: Enter or download patient data
- 2: Check signal quality and perform ECG
- 3: Review ECG before storing, exporting or printing

## THE CARDIOVIT AT-102 G2 DIAGNOSTIC TOOLS

- ⚡ Resting ECG with the SCHILLER ECG analysis program ETM as well as ETM Sport, the automatic interpretation for athletes.
- ⚡ Resting rhythm recording up to 10 minutes.
- ⚡ Optimal signal quality for paediatric ECGs.
- ⚡ Basic exercise ECG (optional) with a large selection of protocols for treadmills and ergometers.
- ⚡ SCHILLER's Culprit Coronary Artery Algorithm (CCAA, optional) locates the obstruction site in the coronary artery and thereby the size of the cardiac area at risk.



The CARDIOVIT AT-102 G2 is also available with a trolley.

## ADD SPIROMETRY TO YOUR CARDIOVIT AT-102 G2

With the SpiroScout SP plus, the CARDIOVIT AT-102 G2 also becomes a spirometer. Thanks to GANSHORN's ultrasound technology, spirometry recordings are highly accurate and the sensor is maintenance- and calibration-free. The disposable ScoutTubes ensure maximum hygiene and patient safety. ScoutTubes are CO<sub>2</sub> neutral.

The spirometry option offers FVC, SVC and MVV with a large selection of reference standards, pre and post tests, as well as interpretation of the FVC results.

The device provides feedback on the procedure's quality during the test and after each trial.

Complete spirometry results can be printed on thermal paper as well as exported to EMR/HIS systems.



An incentive screen is included for trials with children.



**Americas**  
SCHILLER Americas Inc.  
Doral, Florida 33172

**North America:**  
Phone +1 786 845 06 20  
Fax +1 786 845 06 02  
sales@schilleramericas.com  
www.schilleramericas.com

**Latin America & Caribbean**  
Phone +1 305 591 11 21  
Fax +1 786 845 06 02  
sales@schilleramericas.com  
www.schilleramericas.com



**Asia**  
SCHILLER Asia-Pacific / Malaysia  
52200 Kuala Lumpur, Malaysia  
Phone +603 6272 3033  
sales@schiller.com.my  
www.schiller-asia.com



**Austria**  
SCHILLER Handelsgesellschaft m.b.H.  
A-4040 Linz  
Phone +43 732 70 99 0  
Fax +43 732 757 000  
sales@schiller.at  
www.schiller.at



**China**  
Alfred Schiller (Beijing) Medical Equipment  
Co., Ltd.  
100015 Beijing, China  
Phone +86 010 52007020  
Fax +86 010 52007020-8016  
info@schillerchina.com  
www.schillermedical.cn



**Croatia**  
SCHILLER d.o.o.  
10000 Zagreb  
Phone +385 1 309 66 59  
Fax +385 1 309 66 60  
info@schillerzg.hr  
www.schiller.hr



**France**  
SCHILLER Médical  
F-67160 Wissembourg  
Phone +33 3 88 63 36 00  
Fax +33 3 88 63 36 49  
info@schiller.fr  
www.schiller-medical.com



**France (distribution France)**  
SCHILLER France S.A.S.  
F-77608 Bussy St Georges  
Phone +33 1 64 66 50 00  
Fax +33 1 64 66 50 10  
contact@schillerfrance.fr  
www.schiller-france.com



**Germany**  
SCHILLER Medizintechnik GmbH  
D-85622 Feldkirchen b. München  
Phone +49 89 62 99 81 0  
Fax +49 89 609 50 90  
info@schillermed.de  
www.schillermed.de



**India**  
SCHILLER Healthcare India Pvt. Ltd.  
Mumbai - 400 059, India  
Phone +91 22 6152 3333 / 2920 9141  
Fax +91 22 2920 9142  
sales@schillerindia.com  
www.schillerindia.com



**Poland**  
SCHILLER Poland Sp. z o.o.  
PL-02-729 Warszawa  
Phone +48 22 843 20 89 / +48 22 647 35 90  
Fax +48 22 843 20 89  
schiller@schiller.pl  
www.schiller.pl



**Russia & C.I.S.**  
AO SCHILLER.RU  
125040 Moscow, Russia  
Phone +7 (495) 970 11 33  
Fax +7 (495) 956 29 10  
mail@schiller.ru  
www.schiller.ru



**Serbia**  
SCHILLER d.o.o.  
11010 Beograd  
Phone +381 11 39 79 508  
Fax +381 11 39 79 518  
info@schiller.rs  
www.schiller.rs



**Slovenia**  
SCHILLER d.o.o.  
2310 Slovenska Bistrica  
Phone +386 2 843 00 56  
Fax +386 2 843 00 57  
info@schiller.si  
www.schiller.si



**Spain**  
SCHILLER ESPAÑA, S.A.  
ES-28232 Las Rozas/Madrid  
Phone +34 91 713 01 76  
Fax +34 91 355 79 33  
schiller@schiller.es  
www.schiller.es



**Switzerland**  
SCHILLER-Reomed AG  
CH-8912 Obfelden  
Phone +41 44 744 30 00  
Fax +41 44 740 37 10  
sales@schiller-reomed.ch  
www.schiller-reomed.ch



**Turkey**  
SCHILLER TÜRKIYE  
Okmeydani-Sisli - Istanbul  
Phone +90 212 210 8681 (pbx)  
Fax +90 212 210 8684  
info@schiller.com.tr  
www.schiller-turkiye.com



**United Kingdom**  
SCHILLER UK Ltd.  
Coatbridge, ML5 2DA  
Phone +44 1619 072 565  
sales@schilleruk.com  
www.schilleruk.com

Device availability in your market is subject to regulatory approval.

